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1. General points on the medical air station

1.1 Series

The instructions for use apply to medical air stations (RAM).

«RIETTSCHLE MEDICAL AIR»

The yield of the stations below 10 bar is 2 x 7, 2 x 14, 2 x 34, 2 x 44, 2 x 61, 2 x 88 and 2 x 120 m³/hr in 50 Hz.

The compressor specifications are indicated on the sheets 5 and 7.

1.2 Description

The RAM stations consist of 2 compressors (series 5 or 7), a horizontal tank, a programmable automatic control panel, 2 medical air filters and a hygrometer. They can also include a separator/receiver for condensate, optionally with automatic drain and a homogenising tank.

Each compressor is equipped with a non-return valve so that air cannot return once the compressors are stopped. A manual valve is mounted on every device to isolate the compressor for maintenance. The compressor can also be fitted with an automatic valve for constant pressure regulation, on the series 7 this automatic valve is part of the original fittings.

1.2.1 Compressors

They utilise rotary technology with blades and produce pulse-free air and can run for 24 hours a day.

1.2.2 Tanks

Capacity : Horizontal, 500 or 1000 litres, mild steel version, primed and outside painted. Galvanised steel tanks are optionally available.

Every tank has an operating pressure of 10 bar (g) and a proof pressure of 15 bar (g) with safety valve set to 10 bar (g).

Every tank is equipped with a valve of purge and must be used only for pressure.

1.2.3 Control panel

The control panel is designed to manage the whole air station via a programmable controller and is supplied with a display keyboard and control switches. An additional control box manages filtration and accessories.

A cover with main switch protects the panel (see instructions of the board).

1.2.4 Medical air filtration

Every station includes two medical air filtration lines of the MSD type with an adsorption dryer and necessary filters to meet the requirements of norms NFS 90-155 and NFS 90-140.

1.2.5 Device to measure relative humidity

Each station is supplied with a hygrometer to provide continuous measurement of the degree of moisture as defined in the norm NFS 90-140 (alarm signal if dew point reaches - 40°C).

1.2.6 Accessories

Each station can include a homogenising tank to maintain regular oxygen content. A condensate separator for the main tank can be supplied, optionally with automatic condensate drain.

RAM	7.2.500
RAM	14.2.500
RAM	34.2.500
RAM	44.2.1000
RAM	61.2.1000
RAM	88.2.1000
RAM	120.2.1000

BE 60/2

1.7.2001

**Werner Rietschle
GmbH + Co. KG**

Postfach 1260

79642 SCHOPFHEIM
GERMANY

☎ 07622 / 392-0

Fax 07622 / 392300

E-Mail: info@rietschle.com

<http://www.rietschle.com>

Rietschle (UK) Ltd.

Bellingham Way

NEW HYTHE
KENT ME20 6XS
UNITED KINGDOM

☎ 01622 / 716816

Fax 01622 / 715115

E-Mail: info@rietschle.co.uk

<http://www.rietschle.co.uk>

1.3 Application

These air stations are meant for medical air systems. They are designed to work within a range between two set points, minimum and maximum pressure or for continuous running at a constant pressure of 10 bar (g).

 **The air inlet valve must be placed so as to produce the least possible exhaust gas contamination due to internal combustion engines, or any other contamination source.**

The ambient and aspiration temperatures must range from 5 to 25°C. If not, please contact us.

 **The alarm signals must be plugged in and all the safety devices enabled to avoid any risk to staff. A gaseous aid must be installed in case of emergency.**

1.4 Handling and setting up

 **During normal operating, be careful of the hot surfaces (see IS compressors). Sufficient space is necessary on every side of the station to enable maintenance of the different elements.**

The RAM stations can only function correctly in a horizontal position.

The opening and checking of the control panel and board must be performed when switched off and by a skilled worker according to the norm EN 60204.

The board, the switch unit and the filter box must be kept accessible for display, programming and repair.

The medical air filtration and the fittings must be accessible for maintenance.

The air station must be set up in a ventilated room at maximum temperature of 40°C with adequate ventilation. The room must be easily accessible.

 **For correct setting up and running of the station follow the instructions concerning work protection.**

1. Pressure connection at the end of the filtration.
The main air pipe must be of a suitable diameter to avoid any excessive pressure drop.
2. Check the oil level in the compressor (see IS compressor) and add oil if necessary, check each cap is tight.
3. Connect the control board to the power current equipped with a circuit breaker. Check the compatibility to the supply (current, voltage, frequency).

 **The electrical connections and installation should only be performed by a skilled worker according to the norm EN 60204. The main switch for the supply (breaker) must be provided by the user**

1.5 Setting

 **The starting should be realized by Rietschle or a partner of Rietschle.**

Before using check and tighten the board screws.

1. Run the compressors separately and check the direction of rotation: i.e. at the end driving the compressor motor must turn clockwise. If the rotation is correct, the displayed pressure on the gauge immediately rises to 5.5 bar (see IS compressor). Put the devices into operation individually.
2. After a possible change of direction, regulate the required pressure on the casing box (see IS control board).

Factory threshold: 9,5 bar (g)

Constant pressure: $\pm 0,5$ bar

Current pressure range: 8,5 - 9,8 bar (g)

The parameters must not be modified without agreement of Rietschle.

1.6 Maintenance (see operating instruction for compressor, control panel, air filtration and hygrometer)

For servicing the equipment we ask you to use original RIETSCHLE-SPARE PARTS only to guarantee a trouble free operation.

Important: The wearing parts like filter cartridges and oil should be disposed of with due regard to health and safety regulations.

1.7 Appendix

Repair on Site: For all repairs on site an electrician must disconnect the motor so that an accidental start of the unit cannot happen. All engineers must consult the original manufacturer or one of the Partners of Rietschle. The address of Partners of Rietschle can be obtained from the manufacturer on application. After a repair or before re-installation, follow the instructions as shown under the headings "Installations and Initial Operation".

Lifting and Transport: To lift and transport the station the eyebolts on the tank must be used.

Disposal: The wearing parts (as listed in the spare parts lists) should be disposed of with due regard to health and safety regulations.

Modifications: We inform you that any modification of this material is not allowed by the directive 93 / 42 / EEC. It is absolute necessary to consult us for any modification.

2. Operating instruction for compressors

see enclosed documentation

3. Operating instruction for control panels

see enclosed documentation

4. Operating instruction for medical air filtration

see enclosed documentation

5. Operating instruction for hygrometer

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